

# TRADING

## 2000 s.r.l.

### SPECIFICATIONS

#### L CARNITINA BASE

Identification	IR positive test
Appearance	white crystal or crystalline powder
pH	5.5 – 9.5
Specific optical rotation	-29.0° to -32.0°
Sodium	<= 0.10%
Chloride	<= 0.40%
Potassium	<= 0.20%
Heavy metals	<= 10 PPM
Arsenic	<= 1 PPM
Water content	<= 1.00%
Residue on ignition	<= 0.10%
Assay (dried basis )	97.0 – 103.0%
E. Coli	negative
Salmonella	negative
Total plate count	<= 1000 Cfu/g
Yeast & Mold	<= 100 Cfu/g
Conclusion:	Standard of USP32
Storage :	Keep container unopened in cool, dry place, prevent from light

---

Trading 2000 - Società Responsabilità Limitata - Socio Unico

Via Segantini, 69 - 20143 MILANO – Tel. 02.89409101 – Fax 02.8373707

e-mail: [info@trading2000.it](mailto:info@trading2000.it) – <http://www.trading2000.it>

Cod. Fisc. / P.Iva 11636690155 – C.C.I.A.A. Milano 1484779 – Cap. Soc. €102.260,00 int vers. – Reg.Trib. N. 358568